

PET-MPI Annotation Tool

Clinician Study View: Quick Operating Guide

Purpose. This tool records your image-level classification and the grid-cell regions that informed your decision. It is designed for research annotation and does not replace normal clinical interpretation or reporting.

1 Before starting

Use the application in Chrome at **100% browser zoom**. The right panel displays the polar map at its fixed native Study View size of **1024 × 1024**. If the whole map is not visible, use the panel scroll bars; do not change browser zoom during the annotation session. Use the same practical display setup throughout a session.

Confirm that the progress panel on the left shows the expected image number and saved count before beginning work.

2 Annotating one polar map

1. Select a tool	Choose Brush to add grid cells or Eraser to remove selected cells.
2. Choose brush size	Use the compact size list: 1 (8 px) , 2 (16 px) , 4 (32 px) , 8 (64 px) , or 16 (128 px) . Select the smallest size that lets you mark the decision region clearly.
3. Mark the region	Click or drag over the relevant decision-making region. The mark is cell based: each selected cell is either included or not included.
4. Correct if needed	Use Undo , Redo , or Eraser . You may also change the brush colour for easier viewing; the saved scientific mask remains binary.
5. Classify and save	Select Ischemic or Non-ischemic . Use Preview final mask for a clean white-overlay check, then select Save & next image . After the final image, the app confirms that the session is complete.

3 Viewing guides

The **Grid** control shows or hides the annotation-cell boundaries. The small colour square beside it changes the grid colour. The **17-segment** control is an optional visual guide only; it does not constrain what you can select. Its compact selector changes guide strength, and its colour square changes guide colour.

4 Navigation and review

Images are presented in one saved random order for this session. If the app is reopened during the same session, it restores that same order. A new Study Kit/session uses a new random order.

Use **Previous image** to return to an earlier map. Its saved classification and selected cells will be restored. After making any correction, select **Save & next image** again; the existing result for that image is replaced

rather than duplicated. Each save records the image-level class, binary decision-region mask, coloured working overlay, and clean white-overlay preview.

Practical checklist: review the full map before marking, use the smallest practical brush near a boundary, confirm the image-level class, preview when uncertain, and save before leaving an image.