

Medical Image Annotation Tool

Study View: Quick Operating Guide

Purpose. This tool records grid-cell regions on an image and, when required by the study protocol, an optional binary image-level assessment. It is designed for research annotation and does not replace normal clinical interpretation or reporting.

Before starting

Use the application in Chrome at **100% browser zoom**. The annotation workspace displays each input image at its native dimensions. If the full image is not visible, use the workspace scroll bars; do not change browser zoom during an annotation session. Use the same practical display setup throughout a session.

Confirm that the progress panel on the left shows the expected image number and saved count before beginning work.

Annotating one image

Step	What to do
1. Select a tool	Choose Brush to add grid cells or Eraser to remove selected cells.
2. Choose brush size	Use 1 cell , 2 × 2 , 4 × 4 , 8 × 8 , or 16 × 16 . Select the smallest practical size that lets you mark the relevant region clearly.
3. Mark the region	Click or drag over the image region relevant to the study question. The annotation is cell based: each selected cell is either included or not included.
4. Correct if needed	Use Undo , Redo , or Eraser . You may change the brush colour for easier viewing; the saved scientific mask remains binary.
5. Review and save	Use Review output to check the clean white-overlay preview. If the protocol requires an image-level assessment, select Class 1 or Class 0 before saving. Otherwise, leave both options unselected. Select Save & next image . After the final image, the app confirms that the session is complete.

Viewing the grid

The **Grid** control shows or hides the annotation-cell boundaries. The colour square beside it changes the grid colour. The grid is a visual aid and does not change the content of the saved binary mask.

Navigation and review

Images are presented in one saved random order for the session. If the application is reopened during the same session, it restores that same order. A new Study Kit/session uses a new random order.

Use **Previous image** to return to an earlier image. Its saved assessment and selected cells will be restored. After making any correction, select **Save & next image** again; the existing result for that image is replaced rather than duplicated. Each save records the optional image-level assessment, binary region mask, coloured working overlay, and clean white-overlay preview.

Practical checklist: review the full image before marking, use the smallest practical brush near a boundary, review the output when uncertain, and save before leaving an image.