

Medical Image Annotation Tool

A focused local workflow for creating structured, reproducible image-region annotations.

IMAGE DATA STAYS LOCAL**RESEARCH USE SOFTWARE**

The tool saves an analysis-ready binary mask, a visual overlay, a review image, and—only if the study protocol requires it—an optional binary image-level assessment.

AT A GLANCE

Native image dimensions

Each supported image is annotated at its original dimensions. No fixed anatomy, imaging modality, or image size is assumed.

SUPPORTED INPUTS

PNG · JPG · JPEG

01 · PREPARE**Add assigned images**

Place study files in `Study_Data/input_images`.

02 · ANNOTATE**Mark the relevant region**

Use a logical grid, brush, eraser, and review preview.

03 · COLLECT**Retain study outputs**

Collect one independent results folder per session.

**Study-use boundary**

This software supports research annotation. It does not provide a clinical diagnosis, replace clinical reporting, or make treatment recommendations. Use only authorised, de-identified, or otherwise appropriately governed study material.

BEFORE THE FIRST IMAGE

A consistent setup protects the study record

DISPLAY SETUP**Use Chrome at 100% zoom**

Scroll the workspace when needed. Do not resize the source image during the session.

INPUT CHECK**Confirm the assigned files**

Review the expected image count and use only the files approved for this annotation session.

SESSION ISOLATION**One folder, one session**

Keep an independent `Study_Data` folder for each participant or annotation round.

REVIEW HABIT**Save before leaving an image**

Use the final preview when uncertain. Returning to a prior image deliberately replaces its existing output.

SESSION PROGRESS IS PERSISTENT

The application records one randomised image order for the active `Study_Data` folder and restores that order if the session is reopened.

One image, five deliberate actions

Keep the workflow consistent across a session. The session order is randomised once and restored if the application is reopened.

1

Select a tool

Choose Brush to add grid cells or Eraser to remove them.

2

Choose a brush size

Use 1 cell, 2×2 , 4×4 , 8×8 , or 16×16 . Start as small as practical near a boundary.

3

Mark the region

Click or drag over the area relevant to the study question. Each selected cell is included in the binary mask.

4

Correct and review

Undo, Redo, or Erase as needed. Use Review output to inspect the clean white-overlay preview.

5

Save the record

Optionally select Class 1 or Class 0 when defined by the protocol, then choose Save & next image.

OPTIONAL ASSESSMENT**Class 1 / Class 0**

These labels are intentionally study-defined. Leave both unselected unless the protocol requires a binary image-level assessment.

VISUAL CONTROL

Toggle the grid for a clearer view of the source image. The grid helps with reproducible selection; it is not included in the final preview.

WHAT SAVE CREATES**A complete, image-linked output bundle**

binary_data/ · native-dimension 0/1 mask
masked_images/ · transparent annotation overlay
final_preview_images/ · review image
annotation_state/ · selected grid cells
assessments.csv · optional binary assessment

SESSION DISCIPLINE

Use a separate Study_Data folder for each participant or session. To correct a saved image, return with Previous image, edit, and save again; the old output is deliberately replaced rather than duplicated.